

DEDU-8900 (18) Dental Unit



Standard Configuration

4 hole high-speed handpiece tubing with connector	2 set
Low-speed handpiece tubing with connector	1 set
3-way syringe (cold and hot)	2 piece
Adjustable headrest	1 set
New design double armrest	1 set
New design down tray with Alec material with touch	1 set
With autoclavable silicone pad,maximum bearing temperature is 135-degrees;	
6-posistion for handpieces ,syringe, scaler and micro motor	
Water suction and saliva ejector	1 set each
LED operation lamp with sensor(white and yellow light)	1 set
American style lamp arm	1 set
Build-in model LCD X-film viewer	1 set
24V DC motors	1 set
Assistant control system with Alec material with touch	1 set
Water purified system	1 set
Water heating system	1 set
Water bottle	1 pcs

Auto water supply and flush cuspidor system	1 set each
Fixed ceramic cuspidor	1 set
Alec material table	1 set
Multi-function foot pedal DU-106	1 set
Big dentist stool	1 piece
PU leather cushion	1 set
PU Foot Pad	1 set
Solenoid valve DU-097	2 pcs
ABS plastic covers	1 set
PVC air and water tubes	
All kinds of copper valves	
Chair with 9 memories	
With low water alarm function	
Chair with safety system	
Linked movement chair frame	1 set
Aluminium backrest	1 set
New base(steel paste with strong PU cover)	1 set
Packing size: 170*100*120CM	

Options

Handpiece, scaler, curing light, camera system, air compressor, suction unit, oral camera system

DEDU-8900 (18)

Foldable Type Dental Unit



Standard Configuration

4 hole high-speed handpiece tubing with connector	2 set
Low-speed handpiece tubing with connector	1 set
3-way syringe (cold and hot)	2 piece
Strong and weak suction	1 set each
Adjustable headrest	1 set
Double armrest	1 set
Movable tool tray with 9 memories with touch	1 set
With autoclavable silicone pad,maximum bearing temperature is 135°C;	
6-posistion for handpieces ,syringe,scaler and micro motor	
Water suction and saliva ejector	1 set each
Sensor operation sensor lamp DU-152	1 set
24V OPG film viewer	1 set
24V DC motors	1 set

Assistant control system with touch	1 set
Water purified system	1 set
Water heating system	1 set
Water bottle	2 pcs
Auto water supply and flush cuspidor system	1 set each
Turnable ceramic cuspidor	1 set
Foot pedal DU-108	1 set
Big dentist stool	1 piece
PU leather cushion	1 set
Solenoid valve DU-097	2 pcs
Three folding chair frame	1 set
ABS plastic covers	1 set
PVC air and water tubes	
All kinds of copper valves	
Linked movement chair frame	1 set
Chair without safety system	
Packing size: 142*100*120cm, 16unit/20GP, 32unit/40HQ	

Options

Air compressor, suction unit, oral camera system

DEDU-8000 (18) Dental Unit



Standard Configuration

4 hole high-speed handpiece tubing with connector	2 set
Low-speed handpiece tubing with connector	1 set
3-way syringe (cold and hot)	2 piece
Adjustable headrest	1 set
New design double armrest	1 set
Top-mounted instrument tray with 9 memories with touch	1 set
5-position for handpieces, syringe, and scaler	
Stainless steel operation table	1 set
Water suction and saliva ejector	1 set each
LED operation lamp with sensor	1 set
American style lamp arm	1 set
Built in 24V X-film viewer	1 set
24V DC motors	1 set
Assistant control system with touch	1 set
Water purified system	1 set
Water heating system	1 set
Water bottle	1 pcs
Auto water supply and flush cuspidor system	1 set each

Turnable ceramic cuspidor	1 set
Multi-function foot pedal DU-108	1 set
Big dentist stool	1 piece
PU leather cushion	1 set
PU Foot Pad	1 set
Solenoid valve DU-097	2 pcs
ABS plastic covers	1 set
PVC air and water tubes	
All kinds of copper valves	
Linked movement chair frame	1 set
New base(steel paste with strong PU cover)	1 set
Chair without safety system	
Packing size: 1420*1100*1200mm, 16unit/20GP, 32unit/40HQ.	

Options

Air compressor, suction unit, oral camera system

DEDU-9000 (18) Dental Unit



Standard Configuration

4 hole high-speed handpiece tubing with connector	2 set
Low-speed handpiece tubing with connector	1 set
3-way syringe (cold and hot)	2 piece
Adjustable headrest	1 set
Double armrest	1 set
New design down tray with touch	1 set
5-position for handpieces, syringe and scaler	
Water suction and saliva ejector	1 set each
LED lamp with sensor	1 set
24V OPG film viewer	1 set
24V DC motors	1 set
Assistant control system with touch	1 set
Water purified system	1 set
Water heating system	1 set

Water bottle	2 pcs
Tissue box	1 set
Auto water supply and flush cuspidor system	1 set each
Turnable glass cuspidor	1 set
Multi-function foot pedal DU-108	1 set
Dentist stool DU-008	1 piece
PU leather cushion	1 set
Solenoid valve DU-097	2 pcs
ABS plastic covers	1 set
PVC air and water tubes	
All kinds of copper valves	
Without linked movement chair frame	1 set
Chair without memory	
Chair without safety system	
Packing size: 142*100*120cm, 16unit/20GP, 32unit/40HQ	

Options

Handpieces, scaler, curing light, air compressor, Oral camera system

DEDU-2311 (18) Dental Unit



Standard Configuration

4 hole high-speed handpiece tubing with connector	2 set
Low-speed handpiece tubing with connector	1 set
3-way syringe (cold and hot)	2 piece
Adjustable headrest	1 set
Double armrest	1 set
New design down-mounted tool tray with touch	1 set
5-position for handpieces, syringe and scaler	
Water suction and saliva ejector	1 set each
LED reflector operation lamp with sensor	1 set
American style lamp arm	1set
Build in 24V X-film viewer	1 set
24V DC motors	1 set
Asistant control system with touch	1 set
Water purified system	1 set
Water heating system	1 set

Water bottle	1 pcs	
Auto water supply and flush cuspidor system	1 set each	
Turnable ceramic cuspidor	1 set	
Multi-function foot pedal DU-108	1 set	
Dentist stool DU-008	1 piece	
PU leather cushion		1 set
Solenoid valve DU-097	2 pcs	
ABS plastic covers		1 set
PVC air and water tubes		
All kinds of copper valves		
Chair without memory		
Chair frame without linked movement		
Without safety system		
Packing size: 142*100*110cm 16unit/20GP, 32unit/40HQ		

Options

Handpieces, scaler, curing light, air compressor, Oral camera system

DEDU-2305 (18) Dental Unit



Standard Configuration

4 hole high-speed handpiece tubing with connector	2 set
Low-speed handpiece tubing with connector	1 set
3-way syringe (cold and hot)	2 piece
Adjustable headrest	1 set
Double armrest	1 pcs
Down-mounted tool tray	1 set
5-position for handpieces, syringe and scaler	
Water suction and saliva ejector	1 set each
LED lamp with sensor	1 set
24V X-film viewer	1 set
24V DC motors	1 set
Assistant control system with control panel DU-042	1 set
Water purified system	1 set
Water heating system	1 set
Water bottle	1 pcs
Auto water supply and flush cuspidor system	1 set each

Turnable ceramic cuspidor	1 set
Multi-function foot pedal DU-107	1 set
Dentist stool DU-008	1 piece
PU leather cushion	1 set
Solenoid valve DU-097	2 pcs
ABS plastic covers	1 set
PVC air and water tubes	
All kinds of copper valves	
Chair without memory	
Chair without safety system	
Packing size: 142*100*120cm 16unit/20GP, 32unit/40HQ	

Options

Handpieces, scaler, curing light, air compressor, suction unit, oral camera system

DEDU-2305 (17) Dental Unit



Standard Configuration

4 hole high-speed handpiece tubing with connector	2 set
Low-speed handpiece tubing with connector	1 set
3-way syringe (cold and hot)	2 piece
Adjustable headrest	1 set
Double armrest	1 set
New design down-mounted tool tray with touch	1 set
5-position for handpieces, syringe and scaler	
Water suction and saliva ejector	1 set each
LED operation lamp with sensor	1 set
Built in 24V X-film viewer	1 set
24V DC motors	1 set
Assistant control system with touch	1 set
Water purified system	1 set
Water heating system	1 set
Water bottle	1 pcs
Auto water supply and flush cuspidor system	1 set each
Turnable glass cuspidor	1 set

Multi-function foot pedal DU-108	1 set	
Dentist stool DU-008	1 piece	
PU leather cushion		1 set
Solenoid valve DU-097	2 pcs	
ABS plastic covers		1 set
PVC air and water tubes		
All kinds of copper valves		
Chair without memory		
Linked without movement chair frame	1 set	
Without safety system		
Packing size: 142*100*120cm 16unit/20GP, 32unit/40HQ		

Options

Handpieces, scaler, curing light, air compressor, Oral camera system