

IMAR-F50

Immuno Fluorescence



Features

AR-F50 is an advanced on-site analyzer based on immunofluorescence, intending to deliver accurate testing results in minutes merely from 100ul sample. Due to its simple and compact design, FA50 is an ideal point-of-care system for emergency department, intensive care unit.

Product Specifications

Principle: Immunofluorescence quantitative

Display: 5.6 inch touch screen

Sample volume: 100ul

Storage: up to 50,000 results

Calibration mode: Magnetic card

Compact size, easy to carry

Interface: RS232

Reagents

PCT

HS-CRP

cTnl

NT-proBNP

CK-MB

Myo
Cardiac panel (cTnl/CK-MB/Myo)
D-Dimer
HbA1c
T4
TSH
T3*