

ENWG-XP ENT Diagnostic Set



Technical Specifications

1. Ophthalmoscope--- 1.
2. Otoscope Carrier--- 1.
3. Ear Specula--- 4.
4. Nasal Speculum--- 1.
5. Laryngeal Mirror Stem--- 1.
6. Laryngeal Mirrors--- 2.
7. Tongue Depressor Plastic --- 1.
8. Battery Handle--- 1.
9. Spare Bulbs (2 Large and 2. Small 2.5v/0.3A) --- 5.
10. Power Supply: 1.5V X2 "C" Type



ENWG-XPC ENT Diagnostic Set

Technical Specifications

1. Ophthalmoscope--- 1.
2. Otoscope Carrier--- 1.
3. Ear Specula--- 4.
4. Nasal Speculum--- 1.
5. Laryngeal Mirror Stem--- 1.
6. Laryngeal Mirrors--- 2.
7. Tongue Depressor Plastic --- 1.
8. Battery Handle--- 1.
9. Spare Bulbs (2 Large and 2. Small 2.5v/0.3A) --- 5.
10. Transformer--- 1



ENBQZY-I Wall Mounted ENT Diagnostic Set



Features

Wall-mounted, composed of ophthalmoscope, otoscope, nasal speculum, electronic thermometer (Forehead and ear type), storage box and hanging plate, etc. The battery handle is sharing. Light can be adjusted via key control. Power is AC 110~240V, 47~63HZ. Optional parts are tongue depressor, finger clip type pulse oximeter, and blood-glucose meter.

Technical Specifications

| PARAMETER | SPECIFICATION | |
|---|--|----------|
| Ophthalmoscope | | |
| small spot | middle spot | big spot |
| central net | cobalt blue filter | slit |
| Illumination form | small spot, middle spot, big spot, central net, cobalt blue filter, slit | |
| Compensation filter | Polarizing filter, no red green filter | |
| Diopter Compensation | -25~+40 D | |
| Otoscope | | |
| Circle form optical fiber transmission | | |
| Size of the aural speculum | 2.5, 3, 4, 5mm | |
| Magnifier | more than 3 times | |
| Fiber transmission | | |
| Electronic thermometer | | |
| Temperature measurement range | 34°C~43°C | |
| Temperature measurement Precision | ≤ ±0.2°C | |
| Mechanical spring blood pressure meter | | |
| Composed of meter, band and air ball | | |
| Blood pressure gauge: 0~300mmHg | | |
| Upper-arm measurement | manual measurement | |
| Data error | ± 3.75mmHg | |

Storage box

4 kinds of various size aural speculum. Other

consumables can also be stored inside.

ENEJ-II B Otoscope



Technical Specifications

1. Size of the aural: 3, 4, 5, 6, 8mm diameter
2. Power Supply: AC220V;.
3. Illumination Bulb: 2.5V/0.3A.
4. Magnifier: 3 Times.
5. 58x34x22cm.

ENEJ-XPC Otoscope



Technical Specifications

1. Size of the aural: 2.5, 3, 4, 5mm diam
2. Power Supply: AC220V;.
3. Illumination Bulb: 3.5V/0.7A.
4. Magnifier: 3 Times.
5. 36x33x24cm.

Technical Specifications

6. Size of the aural: 2.5, 3, 4, 5mm diam
7. Power Supply: AC220V;.
8. Illumination Bulb: 3.5V/0.7A.
9. Magnifier: 3 Times.
10. 36×33×24cm.

ENEJ-CB OTOSCOPE

Technical Specifications

1. Size of the aural: 2.5, 3, 4, 5mm diam
2. Power Supply: 1.5VX2"AA"Type;
3. Illumination Bulb: 2.5V/0.45A.
4. Magnifier: 3 Times.
5. 39x37x23cm.



ENMHJ-I Fiber Laryngoscope

Technical Specifications

1. Power Supply: 1.5V X2 "C" Type;.
2. Illumination Bulb: 2.5V/0.45A.
3. BIADES: 92, 114, 134mm for adult.



ENMHJ-II Fiber Laryngoscope



Technical Specifications

1. Power Supply: 1.5V X2 "C" Type;.
2. Illumination Bulb: 2.5V/0.45A.
3. BIADES: 75, 100, 125mm for childred.

ENMHJ-III Fiber Laryngoscope



Technical Specifications

1. Power Supply: 1.5V X2 "C" Type,.
2. Illumination Bulb: 2.5V/0.45A.
3. BIADES: 150mm for adult.

EN1800DP

ENT Treating Desk

(for single)



Features

1. Super-mute oil-less host, which greatly reduces the noise and extends the life of the host.
2. The working table which is made of environmental and high-class steel and marble and the baking sheet fission writing desk are smooth, abrasion-resistance, corrosion-resistance, scratch-resistance, impervious and non-stick oil.
3. Spray guns including 6 spray pistols and 2 powder guns, equipped with exquisite loop design and inductive control, can service for a long time.
4. Control panel controls electric chair, film viewer, cold light source and others.
5. Inner Waste tray and wastebin.

Technical Specifications

| PARAMETER | SPECIFICATION |
|---------------------------|-------------------------------|
| Power source | 220AC 50-60Hz |
| Total Power | 1200W |
| Positive Pressure Pump | 100w |
| Negative Pressure Pump | 250w |
| Versatile Cold Light Lamp | 24V 50W Adjustable brightness |
| Table-Board Dimension | L2020×W900×H780 (mm) |

Configuration

| Standard Configuration | quantity |
|--|----------|
| Main Operating Desk | 1 |
| Super-mute Oil-less Positive Pressure Pump | 1 |
| Super-mute Oil-less Negative Pressure Pump | 1 |
| Mute Instrument | 1 |
| Spray Gun | 3 |
| Absorbing Gun | 1 |
| Blowing Gun | 1 |
| Warm-up Demist Device | 1 |
| Versatile Cold Light Lamp | 1 |
| Controller(weak electric controls strong one) | 1 |
| Main and Auxiliary Suction Bottles | 1 |
| Robot Arm and Upright | 1 |
| Equipment Set | 1 |
| Ointment Jar | 4 |
| Wastebin | 1 |
| Suction Tip | 1 |
| Optional | |
| Treating Chair(electric,semi-electric>manual any one is available) | 1 |
| Swivel Chair for Medical | 1 |
| Phial | 4 |

EN1800DH

ENT Treating Desk

(for single)



Features

1. Super-mute oil-less host, which greatly reduces the noise and extends the life of the host.
2. The working table which is made of environmental and high-class steel and marble and the baking sheet fission writing desk are smooth, abrasion-resistance, corrosion-resistance, scratch-resistance, impervious and non-stick oil.
3. Spray guns including 3 spray pistols and 6 powder guns, matched with exquisite loop design, inductive control, long service life.
4. Control panel controls electric chair, film viewer, cold light source and others.
5. Inner Waste tray and wastebin.

Technical Specifications

| PARAMETER | SPECIFICATION |
|---------------------------|-------------------------------|
| Power source | 220AC 50-60Hz |
| Total Power | 900W |
| Positive Pressure Pump | 100w |
| Negative Pressure Pump | 250w |
| Versatile Cold Light Lamp | 24V 50W Adjustable brightness |
| Table-Board Dimension | L1570×W600×H780 (mm) |
| Configuration | |
| Standard Configuration | quantity |
| Main Operating Desk | 1 |

| | |
|--|---|
| Super-mute Oil-less Positive Pressure Pump | 1 |
| Super-mute Oil-less Negative Pressure Pump | 1 |
| Mute Instrument | 1 |
| Spray Gun | 4 |
| Absorbing Gun | 1 |
| Blowing Gun | 1 |
| Warm-up Demist Device | 1 |
| Versatile Cold Light Lamp | 1 |
| Cold Light Source | 1 |
| Controller(weak electric controls strong one) | 1 |
| Main and Auxiliary Suction Bottles | 1 |
| Robot Arm and Upright | 1 |
| Equipment Set | 2 |
| Ointment Jar | 4 |
| Wastebin | 1 |
| Suction Tip | 1 |
| Keyboard Tray | 1 |
| Storage Drawer | 2 |
| Computer Storage Cabinet | 1 |
| Optional | |
| Treating Chair(electric,semi-electric>manual any one is available) | 1 |
| Swivel Chair for Medical | 1 |
| Film Viewer Light | 1 |
| Phial | 3 |
| Endoscope Image Processing Systems(software contained) | 1 |

EN1800SH

ENT Treating Desk

(for Double)



Features

1. Super-mute oil-less host, which greatly reduces the noise and extends the life of the host.
2. The working table which is made of environmental and high-class steel and marble and the baking sheet fission writing desk are smooth, abrasion-resistance, corrosion-resistance, scratch-resistance, impervious and non-stick oil.
3. Spray guns including 3 spray pistols and 6 powder guns, matched with exquisite loop design, inductive control, long service life.
4. Control panel controls electric chair, film viewer, cold light source and others.
5. Inner Waste tray and wastebin.


Technical Specifications

| PARAMETER | SPECIFICATION | |
|---------------------------|-------------------------------|--|
| Power source | 220AC 50-60Hz | |
| Total Power | 1200w | |
| Positive Pressure Pump | 100w | |
| Negative Pressure Pump | 250w | |
| Versatile Cold Light Lamp | 24V 50W Adjustable brightness | |
| Table-Board Dimension | L920×W900×H780(mm) | |
| Configuration | | |
| Standard Configuration | quantity | |
| Main Operating Desk | 1 | |

| | |
|---|---|
| Super-mute Oil-less Positive Pressure Pump | 2 |
| Super-mute Oil-less Negative Pressure Pump | 2 |
| Mute Instrument | 2 |
| Spray Gun(4 straight and 2 bent) | 6 |
| Absorbing Gun | 2 |
| Blowing Gun | 2 |
| Warm-up Demist Device | 2 |
| Versatile Cold Light Lamp | 2 |
| Controller(weak electric controls strong one) | 2 |
| Main and Auxiliary Suction Bottles | 2 |
| Robot Arm and Upright | 2 |
| Equipment Set | 2 |
| Ointment Jar | 6 |
| Wastebin | 2 |
| Suction Tip | 2 |
| Optional | |
| Phial | 6 |

EN1800SC ENT Treating Desk

(for Double)



Features

1. Super-mute oil-less host, which greatly reduces the noise and extends the life of the host.
2. The working table which is made of environmental and high-class steel and marble and the baking sheet fission writing desk are smooth, abrasion-resistance, corrosion-resistance, scratch-resistance, impervious and non-stick oil.
3. Spray guns including 6 spray pistols and 2 powder guns, equipped with exquisite loop design and inductive control, can service for a long time.
4. Control panel controls electric chair, film viewer, cold light source and others.
5. Inner Waste tray and wastebin.

Technical Specifications

| PARAMETER | SPECIFICATION |
|---------------------------|-------------------------------|
| Power source | 220AC 50-60Hz |
| Total Power | 1200w |
| Positive Pressure Pump | 100w |
| Negative Pressure Pump | 250w |
| Versatile Cold Light Lamp | 24V 50W Adjustable brightness |
| Table-Board Dimension | L2020×W900×H780 (mm) |
| Configuration | |
| Standard Configuration | quantity |

| | |
|--|---|
| Main Operating Desk | 1 |
| Super-mute Oil-less Positive Pressure Pump | 2 |
| Super-mute Oil-less Negative Pressure Pump | 2 |
| Mute Instrument | 2 |
| Spray Gun | 8 |
| Absorbing Gun | 2 |
| Blowing Gun | 2 |
| Warm-up Demist Device | 2 |
| Versatile Cold Light Lamp | 2 |
| Controller(weak electric controls strong one) | 2 |
| Main and Auxiliary Suction Bottles | 2 |
| Robot Arm and Upright | 2 |
| Equipment Set | 2 |
| Ointment Jar | 6 |
| Wastebin | 2 |
| Suction Tip | 2 |
| Cold Light Source | 2 |
| Double-bit writing Desk | 1 |
| Computer Storage Cabinet | 2 |
| Keyboard Tray | 2 |
| Optional | |
| Treating Chair(electric,semi-electric>manual any one is available) | 2 |
| Swivel Chair for Medical | 1 |
| Film Viewer Light | 1 |
| Phial | 6 |

ENEN1900B ENT Treating Chair

(manual)

Technical Specifications:

1. Seat height from the ground:460mm.
2. Backrest tilt angle:90° ~ 180° .
3. The rotation of treating chair:360,Can be fixed in any position.
4. Headrest extension:0~200mm.
5. Armrest width can be adjusted in the range of 640mm-680mm.
6. Maximum load:135kg.
7. Size:L640mm×W640mm×H1180mm.



EN1900C ENT Treating Chair

(Semi-Electric)

Technical Specifications:

1. The flat lift range of treating chair:560-760mm(200mm).
2. Backrest tilt angle:90° ~ 180° .
3. The rotation of treating chair:360,Can be fixed in any position.
4. Headrest extension:0~200mm.
5. Armrest width can be adjusted in the range of 640mm-680mm.
6. Maximum load:135kg.
7. Take-off and landing speed:10mm/s
8. Input power: $\leq 100W$
9. Size:L800mmxW640mmxH1180mm
10. The foot switch works in safety extra-low voltage.



EN1900D

ENT Treating Chair

(Electric)



Technical Specifications:

1. The flat lift range of treating chair:560-760mm(200mm).
2. Backrest tilt angle:90° ~ 180° .
3. The rotation of treating chair:360,Can be fixed in any position.
4. Headrest extension:0~200mm.
5. Maximum load:135kg.
6. Take-off and landing speed:10mm/s
7. Input power: $\leq 150W$
8. Size: L790mm×W730mm×H1360mm
9. The control button, and foot switch work in safety extra-low voltage..

EN1900E

Swivel Chair for Medical

(Doctor Chair)



Technical Specifications:

1. Can carry 100kg without deformation and crack phenomenon..
2. Seat lift range: 470mm-600mm.
3. Dimension:L640×W640×H720mm (mm) .
4. With sliding wheels at the bottom, when doctors are operating, it can be moved freely..